



Surgical Uterine Sound (Sterile)

Datasheets





Applications

A surgical uterine sound is a medical instrument used to measure the depth of the uterine cavity. It's a slender, flexible rod with markings indicating depth in centimeters. This tool is commonly used in gynecological procedures like hysteroscopy and endometrial biopsy as listed below:

01. To determine the appropriate length of instruments to be used during procedures like hysteroscopy.
02. To identify any abnormalities or deviations in the uterine structure.
03. To aid in the insertion of other instruments, such as biopsy forceps or intrauterine devices (IUDs).

Specification

Part Number	Length (mm)	Proximal Diameter (mm)	Distal Diameter (mm)	Color
SUS01A1	245	3.5	2.5	White

Packing

Dimensions (mm)			Quantity			Weight / Carton (kg)	
Pouch	Box	Carton	Q'ty / Box	Boxes / Carton	Q'ty / Carton	Net	Gross
60 x 330 x 5	75 x 50 x 331	369 x 352 x 411	25	35	875	3.59	8.17

Instructions for Use

01. Ensure the uterine sound is properly sterilized before use.
02. Gently insert the sound into the cervical canal and advance it into the uterine cavity.
03. Read the depth measurement on the sound to determine the uterine cavity depth.
04. Carefully withdraw the sound from the uterus.

Precautions and Considerations

01. Maintain strict aseptic technique throughout the procedure to prevent infection.
02. Insert the sound slowly and gently to avoid injury to the uterine tissue.
03. Watch for signs of complications, such as bleeding or uterine perforation.
04. Dispose of the used sound according to medical waste disposal guidelines.

It's crucial to note that the use of a uterine sound should be performed by a qualified healthcare professional. Self-use is not recommended and can lead to serious complications.

FDA Listed



LI Medical Corporation

2F, No. 43, Zhongxing Road, Xizhi District,
New Taipei City 221012, Taiwan
Tel : +886-2-8646-1999
Fax : +886-2-2647-1999

Representative / Distributor