



# Pessary – Ring

## Datasheets





## Descriptions

A Ring Pessary is a circular, flexible medical device made of silicone that is inserted into the vagina to support weakened pelvic tissues. It is one of the most commonly used types of pessaries due to its simple design and ease of use. Some versions have a central support membrane for added stability.

## Specifications

| Part Number | Fitting Set | #0   |    | #1   |    | #2   |    | #3   |     | #4   |    | #5   |    | Feature      |
|-------------|-------------|------|----|------|----|------|----|------|-----|------|----|------|----|--------------|
|             |             | Inch | mm | Inch | mm | Inch | mm | Inch | mm  | Inch | mm | Inch | mm |              |
| PR          | None        | 1.75 | 44 | 2.00 | 51 | 2.25 | 57 | 2.50 | 64  | 2.75 | 70 | 3.00 | 76 |              |
| PRS         | Available   | 1.75 | 44 | 2.00 | 51 | 2.25 | 57 | 2.50 | 64  | 2.75 | 70 | 3.00 | 76 | With Support |
| Part Number | Fitting Set | #6   |    | #7   |    | #8   |    | #9   |     | #10  |    |      |    | Feature      |
|             |             | Inch | mm | Inch | mm | Inch | mm | Inch | mm  | Inch | mm | Inch | mm |              |
| PR          | None        | 3.25 | 83 | 3.50 | 89 | 3.75 | 95 | 4.00 | 102 | -    | -  |      |    |              |
| PRS         | Available   | 3.25 | 83 | 3.50 | 89 | 3.75 | 95 | 4.00 | 102 | -    | -  |      |    | With Support |

## Packing

| Part Number | Dimensions (mm) |             |             | Quantity   |              |               | Weight / Carton (kg) |       |
|-------------|-----------------|-------------|-------------|------------|--------------|---------------|----------------------|-------|
|             | Pouch           | Box         | Carton      | Q'ty / Box | Box / Carton | Q'ty / Carton | Net                  | Gross |
| PR00B1      | 120x316x12      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 2.64                 | 8.59  |
| PR01B1      | 120x316x13      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 3.60                 | 9.55  |
| PR02B1      | 120x316x13      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 4.32                 | 10.27 |
| PR03B1      | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 5.04                 | 10.76 |
| PR04B1      | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 5.76                 | 11.48 |
| PR05B1      | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 6.24                 | 11.96 |
| PR06B1      | 120x356x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 6.96                 | 13.23 |
| PR07B1      | 120x356x13      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 3.72                 | 7.78  |
| PR08B1      | 120x356x16      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 6.36                 | 10.42 |
| PR09B1      | 120x356x16      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 6.96                 | 11.02 |
| PRS00B1     | 120x316x12      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 3.36                 | 9.31  |
| PRS01B1     | 120x316x13      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 3.84                 | 9.79  |
| PRS02B1     | 120x316x13      | 180x130x165 | 369x352x411 | 20         | 12           | 240           | 5.04                 | 10.99 |
| PRS03B1     | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 5.52                 | 11.24 |
| PRS04B1     | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 6.24                 | 11.96 |
| PRS05B1     | 120x316x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 6.72                 | 12.44 |
| PRS06B1     | 120x356x13      | 360x130x165 | 369x352x411 | 40         | 6            | 240           | 7.68                 | 13.95 |
| PRS07B1     | 120x356x13      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 4.20                 | 8.26  |
| PRS08B1     | 120x356x16      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 6.48                 | 10.54 |
| PRS09B1     | 120x356x16      | 180x130x165 | 369x352x411 | 10         | 12           | 120           | 7.08                 | 11.14 |

## Applications

Ring pessaries are primarily used to manage:

1. Pelvic Organ Prolapse (POP) : Supports the uterus, bladder, or rectum in cases of mild to moderate prolapse.
2. Stress Urinary Incontinence (SUI) : Helps prevent urine leakage by supporting the urethra.
3. Pregnancy Support : Occasionally used to help prevent preterm labor by providing cervical support.
4. Post-Surgical Support : Can be used after pelvic surgery to maintain structural integrity during healing.

## Precautions

1. Proper Fitting is Essential : A healthcare provider should determine the correct size to prevent discomfort or dislodgement.
2. Regular Cleaning and Maintenance : Patients should be educated on how to remove, clean, and reinsert the pessary, or schedule regular check-ups if unable to manage it themselves.
3. Monitor for Irritation or Infection : Symptoms like vaginal discharge, bleeding, or discomfort may indicate the need for adjustment or removal.
4. Avoid Prolonged Use Without Supervision : Leaving a pessary in for extended periods without cleaning can lead to complications such as infections, ulcers, or tissue damage.
5. Check for Allergic Reactions : Although rare, some individuals may have sensitivities to the materials used in pessaries.

FDA 510(k) K031463



### LI Medical Corporation

2F, No. 43, Zhongxing Road, Xizhi District,  
New Taipei City 221012, Taiwan  
Tel : +886-2-8646-1999  
Fax : +886-2-2647-1999

Representative / Distributor