



Pessary Fitting Sets

Datasheets





Descriptions

Medical-grade silicone **Pessary Fitting Sets** are reusable clinical tools intended to assist healthcare professionals in selecting the appropriate type and size of vaginal pessary for individual patients.

The fitting pessaries are not intended for long-term patient wear and are used exclusively during clinical assessment to determine optimal pessary selection prior to dispensing a patient-use pessary.

Specifications

| Part Number | Description | Size | Packing Dimensions (mm) | Net (kg) | Gross (kg) |
|-------------|--|---------|-------------------------|----------|------------|
| PFSS01C1 | Pessary Fitting Set, Ring without Support | #2 – #7 | 180 x 140 x 110 | 0.149 | 0.989 |
| PFSS02C1 | Pessary Fitting Set, Cube without Drain | #2 – #7 | 200 x 150 x 90 | 0.224 | 1.048 |
| PFSS03C1 | Pessary Fitting Set, Dish with Support | #2 – #7 | 180 x 140 x 110 | 0.170 | 1.010 |
| PFSS04C1 | Pessary Fitting Set, Oval with Support | #2 – #7 | 180 x 140 x 110 | 0.156 | 0.996 |
| PFSS05C1 | Pessary Fitting Set, Gellhorn | #2 – #7 | 270 x 220 x 70 | 0.289 | 1.350 |
| PFSS06C1 | Pessary Fitting Set, Ring with Support | #2 – #7 | 180 x 140 x 110 | 0.165 | 1.005 |
| PFSS07C1 | Pessary Fitting Set, Ring with Knob with Support | #2 – #7 | 180 x 140 x 110 | 0.211 | 1.051 |
| PFSS08C1 | Pessary Fitting Set, MarLand with Support | #2 – #7 | 180 x 140 x 110 | 0.279 | 1.120 |
| PFSS09C1 | Pessary Fitting Set, Dish without Knob without Support | #2 – #7 | 180 x 140 x 110 | 0.145 | 0.986 |

Applications

Pessary Fitting Sets are intended for temporary intravaginal use by qualified healthcare professionals for the purpose of pessary selection.

Clinical Applications Include:

1. Pelvic Organ Prolapse (POP) – anterior, posterior, apical, or combined.
2. Stress Urinary Incontinence (SUI)
3. Mixed Urinary Incontinence (MUI)
4. Symptomatic uterine prolapse or retroversion
5. Postpartum pelvic floor support assessment
6. Pre-surgical or nonsurgical conservative management
7. Patients who are not candidates for surgery or prefer non-invasive treatment

Precautions

1. For use by **trained healthcare professionals only**
2. Fitting pessaries are **not intended for patient take-home use**
3. **Do not use** if the device is managed or shows signs of material degradation
4. **Do not use oil-based lubricants**
5. Use caution in patients with:
 - * Active vaginal infections
 - * Severe vaginal atrophy (unless appropriately managed)
 - * Unexplained vaginal bleeding

FDA 510(k) K031463



LI Medical Corporation

2F, No. 43, Zhongxing Road, Xizhi District,
New Taipei City 221012, Taiwan
Tel : +886-2-8646-1999
Fax : +886-2-2647-1999

Representative / Distributor