



Pessary – Gehrung Datasheets





Descriptions

A Gehrung pessary is a flexible, U-shaped medical device made of medical-grade silicone. It is designed to provide support for more complex pelvic organ prolapse (POP) cases, particularly uterine prolapse with a retroverted (tilted) uterus. Unlike rigid pessaries, the Gehrung pessary can be custom-shaped to fit the contours of the vaginal walls, making it a good option for patients with unique anatomical needs.

Specifications

Part Number	Fitting Set	#0		#1		#2		#3		#4		#5		Feature
		Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	
PGHS	None	2.00	50	2.15	55	2.35	60	2.55	65	2.75	70	2.95	75	With Drain

Part Number	Fitting Set	#6		#7		#8		#9		#10		Feature
		Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	
PGHS	None	3.15	80	3.35	85	3.55	90	3.75	95	-	-	With Drain

Packing

Part Number	Dimensions (mm)			Quantity			Weight / Carton (kg)	
	Pouch	Box	Carton	Q'ty / Box	Box / Carton	Q'ty / Carton	Net	Gross
PGHS00B1	120x316x54	180x130x165	369x352x411	5	12	60	1.50	4.20
PGHS01B1	120x316x54	180x130x165	369x352x411	5	12	60	1.62	4.32
PGHS02B1	120x316x54	180x130x165	369x352x411	5	12	60	1.74	4.44
PGHS03B1	120x316x54	360x130x165	369x352x411	10	6	60	1.92	4.38
PGHS04B1	120x316x54	360x130x165	369x352x411	10	6	60	2.04	4.50
PGHS05B1	120x316x54	360x130x165	369x352x411	10	6	60	2.34	4.80
PGHS06B1	120x356x54	360x130x165	369x352x411	10	6	60	2.46	5.06
PGHS07B1	120x356x54	180x130x165	369x352x411	5	12	60	2.58	5.41
PGHS08B1	120x356x54	180x130x165	369x352x411	5	12	60	2.82	5.65
PGHS09B1	120x356x54	180x130x165	369x352x411	5	12	60	2.94	5.77

Applications

- Moderate to Severe Uterine Prolapse : Ideal for patients with significant uterine descent, particularly those with a retroverted (backward-tilted) uterus.
- Cystocele (Bladder Prolapse) and Rectocele (Rectal Prolapse) : Provides support to the vaginal walls and adjacent pelvic organs.
- Patients Who Need a Customizable Fit : Because it is moldable, it adapts to different vaginal shapes and offers a more personalized fit than other pessaries.
- Non-Surgical Option for POP Management : Suitable for women who cannot or do not wish to undergo surgery.

Precautions

- Requires Professional Fitting : Due to its flexible design, a healthcare provider must adjust and fit the pessary properly to ensure comfort and effectiveness.
- Regular Medical Follow-Ups Needed : It should be checked and cleaned every 4 to 12 weeks to prevent infections and irritation.
- Potential for Vaginal Irritation or Ulcers : Like other pessaries, prolonged use without proper care may cause pressure sores or tissue erosion.
- Not Suitable for Patients with Severe Vaginal Atrophy : Postmenopausal women with thin vaginal tissue may require vaginal estrogen therapy to reduce discomfort.
- Some Patients May Experience Urinary or Bowel Symptoms : If the pessary applies too much pressure on the bladder or rectum, it may lead to difficulty urinating or defecating.
- Difficult for Some Patients to Remove or Clean Independently : Because of its shape and flexibility, some women may need assistance from a healthcare provider for removal and reinsertion.

FDA 510(k) K031463



LI Medical Corporation

2F, No. 43, Zhongxing Road, Xizhi District,
New Taipei City 221012, Taiwan
Tel : +886-2-8646-1999
Fax : +886-2-2647-1999

Representative / Distributor