



Pessary — Shaatz

Datasheets





Descriptions

A Shaatz pessary is a circular, firm yet flexible medical device made from medical-grade silicone. It features a slightly convex shape with a central indentation to aid in placement and retention within the vagina. The design allows it to provide structural support to pelvic organs that have prolapsed.

Specifications

Part Number	Fitting Set	#0		#1		#2		#3		#4		#5		Feature
		Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	
PSH	None	1.50	38	1.75	44	2.00	51	2.25	57	2.50	64	2.75	70	
Part Number	Fitting Set	#6		#7		#8		#9		#10		Feature		
		Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm			
PSH	None	3.00	76	3.25	83	3.50	89	3.75	95	-	-			

Packing

Part Number	Dimensions (mm)			Quantity			Weight / Carton (kg)	
	Pouch	Box	Carton	Q'ty / Box	Box / Carton	Q'ty / Carton	Net	Gross
PSH00B1	120x316x16	180x130x165	369x352x411	10	12	120	0.96	4.74
PSH01B1	120x316x16	180x130x165	369x352x411	10	12	120	1.44	5.22
PSH02B1	120x316x16	180x130x165	369x352x411	10	12	120	2.16	5.94
PSH03B1	120x316x16	180x130x165	369x352x411	10	12	120	3.36	7.14
PSH04B1	120x316x16	180x130x165	369x352x411	10	12	120	4.20	7.98
PSH05B1	120x316x16	180x130x165	369x352x411	10	12	120	4.92	8.70
PSH06B1	120x356x16	180x130x165	369x352x411	10	12	120	6.00	10.06
PSH07B1	120x356x16	180x130x165	369x352x411	10	12	120	6.96	11.02
PSH08B1	120x356x16	180x130x165	369x352x411	10	12	120	7.68	11.74
PSH09B1	120x356x16	180x130x165	369x352x411	10	12	120	9.00	13.06

Applications

The Shaatz pessary is primarily used in the management of:

- Pelvic Organ Prolapse (POP)** : Commonly for **uterine prolapse** and **cystocele** (bladder prolapse), providing internal support to alleviate symptoms.
: Helps reduce urine leakage caused by weakened pelvic floor muscles.
- Stress Urinary Incontinence (SUI)** : Used when surgery is contraindicated or delayed.
- Post-Surgical or Conservative Treatment Option** : In some cases, it is used to manage **cervical incompetence** or to support the uterus in women at risk of preterm labor.
- Support During Pregnancy**

Precautions

- Proper Sizing & Fitting** : A healthcare professional must ensure the correct size for effectiveness and comfort.
- Regular Cleaning & Maintenance** : Patients must be instructed on proper hygiene, including regular removal and cleaning to prevent infections.
: Prolonged or improper use can lead to vaginal irritation, ulcers, infection, or discomfort.
- Monitoring for Complications** : Regular follow-ups are necessary to assess vaginal health and device positioning.
- Follow-up Appointments** : Contraindicated in cases of active vaginal infections, severe vaginal atrophy, or significant ulceration.
- Not Suitable for All Patients** : Water-based lubricants can aid in insertion and removal while preventing tissue irritation.
- Use of Lubricants**

FDA 510(k) K031463



LI Medical Corporation

2F, No. 43, Zhongxing Road, Xizhi District,
New Taipei City 221012, Taiwan
Tel : +886-2-8646-1999
Fax : +886-2-2647-1999

Representative / Distributor